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
Related topics

Angular frequency, characteristic frequency, resonance frequency, torsion pendulum, torsional vibration, torque, restoring torque, damped/undamped free oscillation, forced oscillation, ratio of attenuation/decrement, damping constant, logarithmic decrement, aperiodic case, creeping.

Principle

If an oscillating system is allowed to swing freely it is observed that the decrease of successive maximum amplitudes is highly dependent on the damping. If the oscillating system is stimulated to swing by an external periodic torque, we observe that in the steady state the amplitude is a function of the frequency and the amplitude of the external periodic torque and of the damping. The characteristic frequencies of the free oscillation as well as the resonance curves of the forced oscillation for different damping values are to be determined.

Tasks

A. Free oscillation

1. To determine the oscillating period and the characteristic frequency of the undamped case.
2. To determine the oscillating periods and the corresponding characteristic frequencies for different damping values. Successive, unidirectional maximum amplitudes are to be plotted as a function of time. The corresponding ratios of attenuation, the damping constants and the logarithmic decrements are to be calculated.
3. To realize the aperiodic case and the creeping.

B. Forced oscillation

1. The resonance curves are to be determined and to be represented graphically using the damping values of A.
2. The resonance frequencies are to be determined and are to be compared with the resonance frequency values found beforehand.
3. The phase shifting between the torsion pendulum and the stimulating external torque is to be observed for a small damping value assuming that in one case the stimulating frequency is far below the resonance frequency and in the other case it is far above it.

Theory and evaluation

A. Undamped and damped free oscillation In case of free and damped torsional vibration torques M_1 (spiral spring) and M_2 (eddy current brake) act on the pendulum. We have

$$M_1 = -D^0 \phi \text{ and } M_2 = -C \dot{\phi}$$

ϕ = angle of rotation

$\dot{\phi}$ = angular velocity = $\dot{\phi}$

D^0 = torque per unit angle

C = factor of proportionality depending on the current which supplies the eddy current brake

The resultant torque

$$M = -D^0 \phi - C \dot{\phi}$$

leads us to the following equation of motion:

$$I \ddot{\phi} + C \dot{\phi} + D^0 \phi = 0 \quad (1)$$

I = pendulum's moment of inertia

$\ddot{\phi}$ = angular acceleration

Dividing Eq. (1) by I and using the abbreviations

$$\delta = \frac{C}{2I} \text{ and } \omega_0^2 = \frac{D^0}{I}$$

results in

$$\ddot{\phi} + 2\delta \dot{\phi} + \omega_0^2 \phi = 0 \quad (2)$$

δ is called the "damping constant" and

$$\omega_0 = \sqrt{\frac{D^0}{I}}$$

the characteristic frequency of the undamped system.

The solution of the differential equation (2) is


$$\phi(t) = \phi_0 e^{-\delta t} \cos \omega t \quad (= \phi_0 e^{-\delta t} \cos \omega t) \quad (3)$$

with

$$\omega = \sqrt{\omega_0^2 - \delta^2} \quad (4)$$

Eq. (3) shows that the amplitude $\phi(t)$ of the damped oscillation has decreased to the e -th part of the initial amplitude ϕ_0 after the time $t = 1/\delta$ has elapsed. Moreover, from Eq. (3) it follows that the ratio of two successive amplitudes is constant.

$$\frac{\phi_n}{\phi_{n+1}} = K = e^{\delta T} \quad (5)$$

K is called the "damping ratio" and the quantity

$$\Lambda = \ln K = \delta T = \ln \frac{\phi_n}{\phi_{n+1}} \quad (6)$$

is called the "logarithmic decrement". Eq. (4) has a real solution only if

$$\omega_0^2 \geq \delta^2.$$

For $\omega_0^2 = \delta^2$, the pendulum returns in a minimum of time to its initial position without oscillating (aperiodic case). For $\omega_0^2 < \delta^2$, the pendulum returns asymptotically to its initial position (creeping).

B. Forced oscillation

If the pendulum is acted on by a periodic torque $M_a = M_0 \cos \omega_a t$ Eq. (2) changes into

$$\ddot{\phi} + 2\delta\dot{\phi} + \omega_0^2\phi = F_0 \cos \omega_a t \quad (7)$$

where $F_0 = \frac{M_0}{I}$

In the steady state, the solution of this differential equation is

$$\phi(t) = \phi_a \cos(\omega_a t - \alpha) \quad (8)$$

where

$$\phi_a = \frac{F_0}{\sqrt{\left[1 - \left(\frac{\omega_a}{\omega_0}\right)^2\right]^2 + \left[\frac{2\delta}{\omega_0} \frac{\omega_a}{\omega_0}\right]^2}} \quad (9)$$

and $\phi_0 = \frac{F_0}{\omega_0^2}$

Furthermore:

$$\tan \alpha = \frac{2\delta\omega_a}{\omega_0^2 - \omega_a^2}$$

respectively

$$\alpha = \arctan \frac{2\delta\omega_a}{\omega_0^2 - \omega_a^2} \quad (10)$$

An analysis of Eq. (9) gives evidence of the following:

1. The greater F_0 , the greater ϕ_a
2. For a fixed value F_0 we have:
 $\phi \rightarrow \phi_{\max}$ for $\omega_a = \omega_0$
3. The greater δ , the smaller ϕ_a
4. For $\delta = 0$ we find:
 $\phi_a \rightarrow \infty$ if $\omega_a = \omega_0$

Resonance curves for different dampings.