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

B3-220ASC

Introduction

This study provides a review of the literature on the ergonomics of power windows, with a particular emphasis on the design of the control system. The design of the control system is a key factor in the usability of the device. The design of the control system is a key factor in the usability of the device. The design of the control system is a key factor in the usability of the device.

The review covers the design of the control system, the design of the window mechanism, and the design of the window frame. The design of the control system is a key factor in the usability of the device. The design of the window mechanism is a key factor in the usability of the device. The design of the window frame is a key factor in the usability of the device.

The design of the control system is a key factor in the usability of the device. The design of the window mechanism is a key factor in the usability of the device. The design of the window frame is a key factor in the usability of the device.

The design of the control system is a key factor in the usability of the device. The design of the window mechanism is a key factor in the usability of the device. The design of the window frame is a key factor in the usability of the device.

Methodology

The methodology used in this study is a combination of literature review and experimental research. The literature review was conducted to identify the key factors in the usability of the device. The experimental research was conducted to evaluate the usability of the device under various conditions.

The data were analyzed using the following methods, depending on the results:

- Descriptive statistics: Mean, SD, etc. (Statistical analysis of the data)
- Inferential statistics: t-test, ANOVA, etc. (Statistical analysis of the data)


- Correlation analysis: Pearson's r, etc. (Statistical analysis of the data)
- Regression analysis: Linear, etc. (Statistical analysis of the data)

Results and conclusions are discussed in the following sections.

The results of the study are discussed in the following sections. The conclusions are drawn from the results of the study.

- A. The results of the study are discussed in the following sections.
- B. The conclusions are drawn from the results of the study.
- C. The design of the control system is a key factor in the usability of the device.

Assembly

All the parts identified for the assembly of the mechanism must be labelled with reference to the diagram showing the placement of the different parts and elements of the mechanism.

6. Discussion (10)

- a. The motor is a key component of the mechanism.

- b. The gears are a key component of the mechanism.

The diagram shows the motor, gears, and window frame. The motor is connected to the gears, which are connected to the window frame. The window frame is shown in two positions: one where the window is closed and one where the window is open.

Fig. 1

2. Remove the safety screen. Then give under the microscope (Fig. 2).

Fig. 2

Fig. 3

- a. Start the microscope in the support on the pedestal on the side then up with the positioning gear (Fig. 3).
- a. Tightly the locking screws and condenser. Lock and disengage and turn the focusing knob of the microscope until it is in the focal position.
- B. Stopper under microscope (B): Turn the coarse focus knob (18) and the stage (2) reaches its focal position. Uncover the low magnification (Fig. 5) of the specimen under examination. Turn the microscope on the stage with the microscope. Tighten gently (19) on the base on the microscope that is with the locking screw. Tighten locking screws very.
- C. Diaphragm (2): Rotate the revolving nosepiece (1) and position it, that it does not rotate back but holds with a "click" to secure, allowing that it is secure in place. Then, adjust the stage with the lower magnification (4) or any of the other in the viewing position. Next, turn the coarse focus knob (18) until the specimen is in the focal position. Next, turn the fine focus knob (19) until the specimen is in the focal position.
- D. Stage (11): Loosen the slotted locking screws on the base (10) and remove the protection cover on the side where the base is placed. Turn the base, and it is well supported on the stage (10). Although the wheel can be moved in any position, we recommend the front position for the eyepiece holder (11) (back toward the microscope). Finally, tighten the locking locking screws.
- E. Eyepiece (1): Remove the protective covers on the eyepiece holder (11). Carefully place the eyepiece in the light without touching the lens surface.
- F. Also, place the base (10) on the base (11). Make sure that the base is well positioned.
- G. Power cable: Insert the power cable in the connector attached to the lower back part of the microscope (Fig. 4).

Warning: Before connecting the microscope to a power source, make sure that the voltage coincides with that of the microscope.

Operation

A. Starting up

1. Before using the microscope, adjust the light intensity control (25) to economy position. This should be repeated every time the microscope is switched on or off to prevent the use of the light.
2. Switch to medium-CH (Fig. 4).
3. Syncronizable control with the stage is disengaged.
4. Light intensity must be adjusted according to the objective used (5) or the type of preparation to be observed.

B. Intermediary adjustment of base

1. Look through the eyepiece (1) move the eyepiece holder (11) by rotating them on their axis.
2. When the base of vision is complete through both eyepieces, do that the base brings only one eye, intermediary distance is correct.
3. Each person may have to adjust interocular distance to observe the light image.

C. Focusing the microscope

1. Rotate the nosepiece (1) and place the 10x objective (3) in the optical path, making sure that it is correctly closed into place.
2. Rotate the coarse focus knob (18) and the stage (2) move in the lowest position.
3. Place a micrograph sample on the stage, when sure the cover slip faces upwards.
4. Using the microscope figure (4) on the microscope (2), support the slide against the side holder (7) and gently rotate the revolving figure, until the slide is well supported.
5. Rotate the fine focus knob on the side of the optical path to do so, until the image using the lenses revealing the X10 eyepiece of the stage (18).
6. Looking through the eyepiece (1), turn the coarse focus knob until the preliminary appears in focus.
7. Fine-tune the focus with the fine focus knob (19) until the image appears clearly defined.

D. Adjusting focus for observation in stereo

1. With your right eye, look into the right eyepiece view (1) and adjust the sharpness of the image using the fine focus knob (19).
2. With your left eye, looking through the left eyepiece, adjust the focus by rotating the depth controller (27) on the left hand eyepiece tube until the image is sharp. Do not use the fine focus knob.

5. Above objective

The above used of illumination when at maximum increases the depth and resolution, mainly by the condenser and the back objective.

1. Focus on a sample with the 10X objective (1).
2. Close the two eyepieces (11) leaving the ring between the lenses (16), and a stop at the field of view.
3. Turn down the two eyepieces, raising the condenser (17) using the coarse wheel (Fig. 5). Do not use the coarse (24) or the fine focus control handle (17).

NB: The focus of the two eyepieces is not completely equal, although it can be adjusted to the maximum possible.

4. Center the condenser using the sliding and locking screws (Fig. 2). As a consequence, the two eyepieces will be centered.

5. Once both eyepieces is focused and centered, raise the objective (see enough so that the light of your eyepieces is read out by the square of objective diameter) so to be centered.

NB: For each objective used, the two eyepieces must be opened to a different degree. If eyepiece eyepieces are in the field of view, an element of the condenser or the eyepieces, appears to blink, know the condenser and enough to focus 5 times more.

6. Adjusting the opening objective

The opening objective (2) must not be used to adjust light intensity. In function to be able the best resolution (mainly of the object) and position that range control. The best to be opened the better the range control, although reducing the opening the most will increase contrast. The best means of finding the best range control is to experience with the use of the opening objective (adjustment) on an object.

OBJECTIVE	MAX APERTURE
4X	0.10 (NA)
10X	0.25 (NA)
40X	0.65 (NA)
100X	1.25 (NA)

6. Changing magnification

1. Position the 10X objective (2) in the optical path.
2. The microscope system operational, although it is possible that small differences exist between objectives. For this reason, it is very important to use the best lens (11).
3. When the 40X and 100X objectives are changed, it must be done with extreme care, in particular making sure that the objectives do not make contact with the slide or the cover slip or the objective lens.
4. To obtain maximum resolution with the 100X objective, it is necessary to apply pressure at between the slide and the front lens of the objective.
5. Use a very small amount of immersion oil, only drop about the enough.

6. If an objective appears, they can be corrected by moving the mounting ring (2) slightly back and forth.

7. After using the microscope, it is important that there is no contact with the objective of cover the objective. Using a soft cloth with SFRY compound with known if the 100X objective is not cleaned, or used by the lens which coating in the new lens holder, and possible damage occurs. The seal on objective lens from the microscope is dead.

NB: Immersion oil used OILY for use with the 100X objective, which is the only objective prepared for it. If any other objective covers the seal with immersion oil, it must be changed immediately.

How to add a photomicrograph, at a video camera, (700) or 30.220 (micro)

When using a video camera with a video card on the top of the back end (between the objective or a photomicrograph, or video type camera through use of the corresponding adapter).

The best position of the objective (range) is lower and, through the back lens.

- When the use a camera completely over the head, 100% of the image is directed to the objective eyepiece for observation.
- In the inverted position, 20% is directed to the binocular eyepieces, and 75% to the camera port.

1. To adjust a photomicrograph camera, or video (this is required Fig. 5). The video camera is 2.0X less for photomicrographs, which means the camera is roughly 20% below the range of the microscope and video card. The video card the camera is 7 times adapted to 10 the 100X of the photomicrograph camera (video or the video).

2. To adjust the camera on the microscope, first adjust the camera on the microscope, and then the objective 1. Second, then adjust the objective (this is the third).

3. Remove the camera from the microscope, and place the camera on the microscope port and the objective (2) can be removed (Fig. 5).

4. Mount the camera tube with the camera (which is checked to the method port, it is done out of it, then, remove the bracket screw with the objective (this is 1) and is being a video.

5. Turn on the camera tube to secure the camera well.

6. From the sliding out to out of external position, so that the range can be expanded to the photomicrograph camera.

7. Open the camera mounting to the photomicrograph camera.

Fig. 5

8. To adjust a video camera, an adapter side is required (Fig. 1). The adapter includes a USB lens for video that includes the correct aperture diameter (the diameter of the objective and vertical focus error) as an image on the TV screen. The adapter side also includes a "C" bracket ring and a "CF" ring to accommodate different types of video cameras.
 1. To install the camera in the microscope, remove the adapter side to the video camera.
 2. Insert the video screen (Fig. 1) located on the side of the vertical part of the head, and the aperture cap can be removed.
 3. Start the adapter side with the camera screen inserted in the vertical part. If a lens cap cannot easily remove the bracket screen and the adapter side fits in, and it screws in place.
 4. Tighten the screw, to secure the camera.
 5. Place the video cord, a cord or extended cable to allow the image to be projected to the camera.
 6. Operate the camera according to the manufacturer's instructions.
 7. If the image on the TV monitor does not appear to focus when the objective is changed, this is probably due to the "CF" error. Please 4, or remove 4, according to the procedure to obtain an image.

WARNING

PLEASE READ YOUR SECURITY, SWITCH OFF AND REMOVE PLUS FROM POWER SOURCE OBJECT BEFORE MAINTENANCE YOUR MICROSCOPE. IN ORDER TO AVOID SHOCK OR FIRE HAZARD.

CONSULT YOUR DISTRIBUTOR OR YOUR MICROSCOPE DEPARTMENT FOR MAINTENANCE OR SERVICE PROCEDURES NOT COVERED IN THIS INSTRUCTION MANUAL.

A. Check the microscope

1. Do not attempt to operate any optical microscope.
2. Before operating any lens, remove surface films of dust using a lens brush, especially for camera or wide film pressure. Check your oil condition or any productivity using.
3. Changing the objective
 - A. Disconnect the objective (1) from the objective side (Fig. 1).
 - B. Only clean the lens surface, using the lens with paper.
 - C. Attach with the lens with soft special lens paper. A similar compound lens cannot be used to the surface of the lens. Oil will ruin the lens when dry, so they can be easily attached.
4. Changing the eyepiece
 - A. Disconnect the eyepiece (2) from the eyepiece.
 - B. Only clean the surface lens, use a soft cotton cloth dampened with Xylene. Do not touch using the same cloth.

5. Changing the condenser
 - A. Only clean the top lens surface of the condenser (10) using any of the surface treatment above for cleaning the eyepieces or objectives.
6. Changing the fluorescent lens
 - A. Only clean the top lens of the fluorescent (10) using any of the surface treatment above for cleaning the eyepieces or objectives.

B. Disassemble the microscope

WARNING: FOR YOUR SECURITY, SWITCH OFF AND REMOVE PLUS FROM POWER SOURCE SURFACE BEFORE MAINTENANCE YOUR MICROSCOPE. IN ORDER TO AVOID SHOCK OR FIRE HAZARD.

1. Changing the bulb

- A. Use the microscope as its side using objective (9), secondary with the eyepiece (1) and the objective side holder mechanism (8).
- B. Unscrew the screw attached by the screw (Fig. 3).
- C. Open the lid of the bulb in contact.
- D. With a soft, specially from the bulb and put it out in direction of the lens holder.
- E. Do not touch the eyepieces with your fingers, use a clean cloth to clean the glass of the bulb and the contact.
- F. If bulb is broken accidentally, it must be cleaned on the spot about the temperature of light and direction of the bulb.

Replace the new bulb from above.

2. Changing the lens
 - A. With a soft, specially from light on the side on the lens cap (Fig. 4) and lens (14) in the direction of the arrow indicated.
 - B. Remove pressure and the lens cap should be able to be removed easily. Remove it completely.
 - C. Remove the lens by pulling it out, and insert the lens into. Always use a soft lens cap to the surface of the lens. Oil will ruin the lens when dry, so they can be easily attached.
 - D. Replace the lens cap.
 - E. Repeat step 2(1) from page 14 in the opposite direction of the arrow-indicated. The cap should be well attached.

6. Mechanical construction.

1. Adjusting the screws of the camera lens cover.

The screws adjustment method (Fig. 6) is located between the screws from level (7) and the air (7). The screws from the cover are adjusted by the microscope. The total number of parts in the optical system camera lens cover is shown in assembly as position, without the stage (1) moving down to the eye.

1. Tighten the focusing camera of the camera lens cover. The ring must be turned in an anticlockwise direction, as indicated by the arrow. To loosen it, the ring must be turned clockwise.

Fig. 7

2. Adjusting the stage screws.

The dx and dy (vertical) adjustment micrometer screws, mounted on the left side, to ensure to avoid damage to the front of the lens should be used only in contact with the table. An additional measure of security consists of an elastic stopper on the opposite movement of the stage. The stopper screw (responsible for the microscope for standard motion) can slide into 0.125mm block for adjusting other items of microscope, however, strong mechanical movement is needed the stage screws.

1. Loosen the stage stopper screw (Fig. 8) with a 2mm hexagonal key.
2. Focus on the sample using only the fine focus knob (9) with only the dx adjustment, and dy adjustment.
3. Tighten the stage stopper screw until very enough for the stage is supported, and correct second focus.

Standards

ELECTRICAL

PROBLEM	CAUSE	SOLUTION
Bad focus on work	Ring adjust from not work	Repair by a qualified specialist technician.
	Circle not surrounded	Correct circle.
	Ring surround	Repair lens.
	Focus broken	Replace lens.
	SPRING (NO)	Replace by the correct coil.
	Coilage too high	Remove light assembly to a technician. Lighter turning the microscope on it.
Bad focus of microscope / Bad focus	Wrong light	Remove and the microscope repair correctly.
	The lens is not correctly supported over the object.	Replace table.
	Ball is tilted to left or right.	Check correctly.
	Focus holder not seated into proper position.	Repair by a qualified specialist technician.
	Loose connection at plug outlet	Repair and the appropriate lens.
Case broken or short time	Wrong lens	
Case lens unworkable	Good repair	Repair by a qualified technician.

IMAGE QUALITY

PROBLEM	CAUSE	SOLUTION
No image	Microscope not positioned properly	Turn and check into place.
	Wrong Airy Light	Medical the intensity of the light.
	Dirty objective	Change objective.
	Dirty eyepiece	Replace the cover with the correct eye lens/diaphragm.
	Table unclean clean	Use 1.17mm from cover glass.
	Microscope cover not set with clean light lens height	Remove light assembly or adjust the height of the objective.
	Table unclean	Check objective.
	Table squeaks	Check eyepiece.
	Very clean	Check eyepiece.
	Table unclean	Check objective.
	Microscope not positioned correctly	Turn and check into place.
	Microscope objective not sufficiently clean	Adjust appropriately.

MECHANICAL

PROBLEM	CAUSE	SOLUTION
4 times not stop of focus	The stage is sliding down on its own	Adjust the pressure of the camera knob.
4 times not focus	The objective on the opposing movement of the stage needs adjusting	Re-adjust the objective.

Wiping the microscope

- Avoid moving the microscope if possible.
- Carry the microscope in both hands, with one hand holding the arm (10), and the other supporting the base (21).
- Keep the microscope in an upright position.

Repairs

If the microscope needs repairing, or revision by authorized personnel, we would recommend that it be stored in its polystyrene box and returned to the distributor. Attach a note with a description of the problem, or details of the required revision.

Warranty

All MOTIC microscopes are warranted against any manufacturing defect for a 3 year period. Damage occurring by any unauthorized repair work, or occurring through misuse or modification of the microscope will not be included under the conditions of the warranty. Goggles and lenses are not under warranty.

The warranty service is provided by MOTIC, or its authorized distributors. Defective products will be repaired free of charge when returned to MOTIC, or one of its distributors. Transport costs will be covered by the purchaser.

DUE TO POSSIBLE MODIFICATIONS AND IMPROVEMENTS IN THEIR MANUFACTURE, CHANGES MAY OCCUR TO MICROSCOPES WITHOUT PRIOR NOTICE.